

'Trust Shipping' is the practice name for 'Trusted Shipping Limited'

Trust Shipping

Hong Kong

T: +852 2910 7722

E: service@trustshipping.org

U: www.trustshipping.org



Breath Alco Analyser

Trust Shipping, as agents, is proud to introduce you, to our 'Breath Alco Analyser'

Our Breath Alco Analyser, are a cost effective method of monitoring for alcohol consumption in your Zero Tolerance testing program.

Our Breath Alco Analyser meet the requirements of testing by the US Department of Transportation (DOT) & US Food & Drug Administration (FDA).



Specifications

- Dimensions 4.75 X 1.75 X / .75 inches (12.1 cm x 4.4 cm x 1.9 cm)
- Weight: 1.7 oz (48 g) without batteries;
- Housing: Shock resistant, molded plastic
- Power: 2 AA Batteries
- Battery Life: 150 Tests
- Sensor: Sensitive Semi-conductor Sensor
- Blowing time: 5 seconds
- Response time: 5 seconds
- Digital display: 2 Digits (.00% B.A.C)
- Operating range:10-40c (50-104f)
- Detection Range: .00% - .40% B.A.C
- Air Sample: 5 seconds
- Operation: 1 Button
- Limited Warranty: 1 Year
- Accuracy: +/- 0.01% at 0.02% B.A.C

Features

- Simple & easy to use
- Light weight & compact
- Folding mouthpiece design
- Two digit LED display
- Beeps when ready to use
- Beeps after breath sample obtained
- Semi-conductor oxide sensor
- Electric pressure sensor to improve breath capture quality