



Trust Shipping

Hong Kong

T: +852 2910 7722

E: service@trustshipping.org

U: www.trustshipping.org

Drug Testing Kits

Trust Shipping, as agents, is proud to introduce you, to our rapid drug testing kits which are the best drug tests on the market today! Quality you can trust at very affordable prices. The test method is simple & self-contained. The Drug Cup results are visible within minutes. It is simple to use, with a positive result indicated with one line and a negative result with two lines.



Our drug test cup method is an innovative drug test cup design that eliminates urine handling and donor tampering, providing a unique and efficient approach for on-site testing. Because of its compact design, round shape this drug test cup is easy to use for men and women.

Our tests are the same tests used by hospitals and clinics. Our test cup method is FDA (US Department of Health & Human Services ~ US Food & Drug Administration) cleared, easy to use and easy to read. However, the panel of drugs which can be tested, have been tailor made for Trusted Shipping Limited, in accordance with the drug panel recommended by the OCIMF and as required by the US DoT.

It is our commitment to ensure our clients go one step beyond mandatory compliance requirements of the US CFR & recommendations of the OCIMF.

Mandatory requirements:

46 CFR 4.06-5: Requirements for alcohol and drug testing following a serious marine incident (effective June 20, 2006). In this regard drug testing must be conducted on each individual engaged or employed on board the vessel who is directly involved in the SMI (Serious Marine Incident).

(i) The collection of drug-test specimens of each individual must be conducted within 32 hours of when the SMI occurred, unless precluded by safety concerns directly related to the incident.

(ii) If safety concerns directly related to the SMI prevent the collection of drug-test specimens from being conducted within 32 hours of the occurrence of the incident, then the collection of drug-test specimens must be conducted as soon as the safety concerns are addressed.

Instructions for testing of any combination of the following drugs:

**AMP1000 COC300 THC50 OPI2000 MET1000 OXY100 BZO300 PCP25 BAR300 MTD300 BUP10
MDMA1000 - Using One Step DOA Cup (Urine)**

AA: INTENDED USE & SUMMARY

The Rapid Response™ Multi-Drug Urine Test Cup is a rapid, qualitative test for the simultaneous detection of multiple drugs and metabolites in human urine at the cut off concentration of

Test	Target Drug	Cut-off
Amphetamine (AMP)	d-Amphetamine	1000 ng/mL
Cocaine (COC)	Benzoyllecgonine	300 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9-COOH	50 ng/mL
Opiates (OPI)	Morphine	2000 ng/mL
Methamphetamines (MET)	Methamphetamine	1000 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL
Barbiturates (BAR)	Secobarbital	300 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Buprenorphine (BUP)	BUO-3-D-Glucuronide	10 ng/mL
Ecstasy (MDMA)	Ecstasy	1000 ng/mL