



iRupt Jr.

Mini Cell Homogenizer




Max. 4000 RPM, 3 x 2 ml

- ♦ Variable speed from 2800 to 4000 RPM
- ♦ Pulse mode function
- ♦ Variable cyclic intervals
- ♦ Compact and light weight
- ♦ Powerful motor to homogenize most challenging samples
- ♦ Lid lock safety

iRupt 24P

Cell / Tissue Homogenizer



Max. 4550 RPM

- ♦ Brushless motor
- ♦ Microprocessor controlled
- ♦ Variable speed from 490 to 4550 RPM
- ♦ Save up to 99 programs
- ♦ Variable cyclic intervals
- ♦ Accommodates different tubes 24 X 1.5 / 2 ml & 12 X 7 ml